



dnaN Antibody

Product Code	CSB-PA326890XA01SYN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29438
Immunogen	Recombinant Serratia marcescens dnaN protein
Raised In	Rabbit
Species Reactivity	Serratia marcescens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Serratia marcescens
Lead Time	Made-to-order (14-16 weeks)
Target Names	dnaN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.