



dsbA Antibody

Product Code	CSB-PA671529XA01DPX
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q44504
Immunogen	Recombinant Azotobacter vinelandii dsbA protein
Raised In	Rabbit
Species Reactivity	Azotobacter vinelandii
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Azotobacter vinelandii
Lead Time	Made-to-order (14-16 weeks)
Target Names	dsbA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.