



Uncharacterized protein gp14 Antibody

Product Code	CSB-PA329504XA01FSP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P24795
Immunogen	Recombinant Escherichia phage D108 (Bacteriophage D108) Uncharacterized protein gp14 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage D108 (Bacteriophage D108)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage D108 (Bacteriophage D108)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.