



sdC Antibody

Product Code	CSB-PA326287XA01SMK
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P26295
Immunogen	Recombinant Streptococcus dysgalactiae subsp. equisimilis (Streptococcus equisimilis) sdc protein
Raised In	Rabbit
Species Reactivity	Streptococcus dysgalactiae subsp. equisimilis (Streptococcus equisimilis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Streptococcus dysgalactiae subsp. equisimilis (Streptococcus equisimilis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	sdc
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.