



ea8.5 Antibody

Product Code	CSB-PA356137XA01FSF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P03755
Immunogen	Recombinant Escherichia phage lambda (Bacteriophage lambda) ea8.5 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage lambda (Bacteriophage lambda)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage lambda (Bacteriophage lambda)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ea8.5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.