



74 Antibody

Product Code	CSB-PA954408XA01FRY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B0FIN6
Immunogen	Recombinant Escherichia phage Phieco32 (Escherichia coli phage phi32) 74 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage Phieco32 (Escherichia coli phage phi32)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage Phieco32 (Escherichia coli phage phi32)
Lead Time	Made-to-order (14-16 weeks)
Target Names	74
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.