



estB Antibody

Product Code	CSB-PA690460XA01FFT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q53547
Immunogen	Recombinant Pseudomonas fluorescens estB protein
Raised In	Rabbit
Species Reactivity	Pseudomonas fluorescens
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pseudomonas fluorescens
Lead Time	Made-to-order (14-16 weeks)
Target Names	estB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.