




fepA Antibody

Product Code	CSB-PA361633XA11ENV
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P05825
Immunogen	Recombinant Escherichia coli (strain K12) fepA protein (23-746aa)
Raised In	Rabbit
Species Reactivity	Escherichia coli (strain K12)
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:2000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia coli (strain K12)
Target Names	fepA

Image



Western Blot

Positive WB detected in: CSB-PA361633XA11ENV(40ng), CSB-PA361633XA11ENV(20ng), CSB-PA361633XA11ENV(10ng), CSB-PA361633XA11ENV(5ng)

All lanes: fepA antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/40000 dilution

Predicted band size: 87 kDa

Observed band size: 100 kDa

Exposure time: 120s

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.