



fla Antibody

Product Code	CSB-PA972141XA01FTR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53606
Immunogen	Recombinant Tatlockia micdadei (Legionella micdadei) fla protein
Raised In	Rabbit
Species Reactivity	Tatlockia micdadei (Legionella micdadei)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Tatlockia micdadei (Legionella micdadei)
Lead Time	Made-to-order (14-16 weeks)
Target Names	fla
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.