



Recombinant Cyprinus carpio Proto-oncogene c-Fos (fos)

Product Code	CSB-EP008790EQE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P79702
Product Type	Recombinant Protein
Immunogen Species	Cyprinus carpio (Common carp)
Purity	>85% (SDS-PAGE)
Sequence	MMFTSLNADC DASSRCSTAS AAAESVACYP LNQTQKFTEL SVSSASFVPT VTAISSCPDL QWMVQPMVSS VAPSNGGARS YNPNPYPKMR VTGTKSPNSN KRARAEQLSP EEEEEKRVRR ERNKMAAAKC RNRRRELTDT LQAETDELED EKSALQNDIA NLLKEKERLE FILAAHKPIC KIPSSSVSPI PAASVPEIHS ITTSVVSTAN APVTTSSSSS LFSSTASTDS FGSTVEISDL EPTLEESLEL LAKAELETAR SVPDVLSSS LYARDWESLY TPANNDLEPL CTPVVTRTPA CTTYTSSFTF TYPENDVFPS CGPVHRRGSS SNDQSSDSL SPTLLTL
Source	E.coli
Target Names	fos
Protein Names	Recommended name: Proto-oncogene c-Fos Alternative name(s): Cellular oncogene fos
Expression Region	1-347
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.