



# fixA Antibody

<b>Product Code</b>	CSB-PA487727XA01EOK
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	B7MAG4
<b>Immunogen</b>	A synthesized peptide derived from Escherichia coli O45:K1 (strain S88 / ExPEC) fixA protein
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Tested Applications</b>	ELISA, Dot blot (guaranteed specific to the immunizing peptide)
<b>Form</b>	Liquid
<b>Conjugate</b>	Non-conjugated
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	Antigen Affinity Purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Product Type</b>	Polyclonal Antibody
<b>Immunogen Species</b>	Escherichia coli O45:K1 (strain S88 / ExPEC)
<b>Lead Time</b>	Made-to-order (14-16 weeks)
<b>Target Names</b>	fixA
<b>Product Modify</b>	<ul style="list-style-type: none"> <li>①. 100ug * KLH-conjugated peptide paired with 100ug * BSA-conjugated peptide (positive control);</li> <li>②. 1ml*pre-immune serum (used as negative control) for free;</li> <li>③. Rabbit polyclonal antibodies purified by Antigen Affinity;</li> </ul>
<b>Product Precision</b>	<ul style="list-style-type: none"> <li>1. CUSABIO guarantees the antibody purity &gt; 95% confirmed by SDS-PAGE.</li> <li>2. CUSABIO guarantees the antibody titer &gt; 1: 64,000 confirmed by ELISA.</li> <li>3. CUSABIO guarantees to offer Peptide-specific Dot blot.</li> <li>4. CUSABIO guarantees 100% risk-free: No qualified antibody, no charge.</li> </ul>
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.