



ftsZ1 Antibody

Product Code	CSB-PA326965XA01RKU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30327
Immunogen	Recombinant Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti) ftsZ1 protein
Raised In	Rabbit
Species Reactivity	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhizobium meliloti (strain 1021) (Ensifer meliloti) (Sinorhizobium meliloti)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ftsZ1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.