



III Antibody

Product Code	CSB-PA526801XA01ELA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O80297
Immunogen	Recombinant Enterobacteria phage If1 (Bacteriophage If1) III protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage If1 (Bacteriophage If1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage If1 (Bacteriophage If1)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	III
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.