



GP90 Antibody

Product Code	CSB-PA343437XA01FTS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P49714
Immunogen	Recombinant Escherichia phage K1E (Bacteriophage K1E) GP90 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage K1E (Bacteriophage K1E)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage K1E (Bacteriophage K1E)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	GP90