



fre Antibody

Product Code	CSB-PA116088XA01VFD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P43127
Immunogen	Recombinant Vibrio harveyi (Beneckea harveyi) fre protein
Raised In	Rabbit
Species Reactivity	Vibrio harveyi (Beneckea harveyi)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vibrio harveyi (Beneckea harveyi)
Lead Time	Made-to-order (14-16 weeks)
Target Names	fre
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.