



ftsZ Antibody

Product Code	CSB-PA323812XA01BRJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P17865
Immunogen	Recombinant Bacillus subtilis (strain 168) ftsZ protein
Raised In	Rabbit
Species Reactivity	Bacillus subtilis (strain 168)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Bacillus subtilis (strain 168)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ftsZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.