



fuml Antibody

Product Code	CSB-PA184767XA01FTT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	D2D3B2
Immunogen	Recombinant Sphingopyxis macrogoltabida (Sphingomonas macrogoltabidus) fuml protein
Raised In	Rabbit
Species Reactivity	Sphingopyxis macrogoltabida (Sphingomonas macrogoltabidus)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Sphingopyxis macrogoltabida (Sphingomonas macrogoltabidus)
Lead Time	Made-to-order (14-16 weeks)
Target Names	fuml
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.