



VII Antibody

Product Code	CSB-PA300852XA01ECJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P69533
Immunogen	Recombinant Enterobacteria phage fd (Bacteriophage fd) VII protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage fd (Bacteriophage fd)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage fd (Bacteriophage fd)
Lead Time	Made-to-order (14-16 weeks)
Target Names	VII
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.