



ftsZ Antibody

Product Code	CSB-PA860156XA01WAN
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ALA3
Immunogen	Recombinant Wigglesworthia glossinidia brevipalpis ftsZ protein
Raised In	Rabbit
Species Reactivity	Wigglesworthia glossinidia brevipalpis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Wigglesworthia glossinidia brevipalpis
Lead Time	Made-to-order (14-16 weeks)
Target Names	ftsZ
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.