



GAPN Antibody

Product Code	CSB-PA305504XA01EWE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P81406
Immunogen	Recombinant Pisum sativum (Garden pea) GAPN protein
Raised In	Rabbit
Species Reactivity	Pisum sativum (Garden pea)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pisum sativum (Garden pea)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GAPN
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.