



37 Antibody

Product Code	CSB-PA863874XA01FSC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9G0B5
Immunogen	Recombinant Escherichia phage AR1 (Bacteriophage AR1) 37 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage AR1 (Bacteriophage AR1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage AR1 (Bacteriophage AR1)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	37