



fixB Antibody

Product Code	CSB-PA345693XA01CME
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P53578
Immunogen	Recombinant Clostridium saccharobutylicum fixB protein
Raised In	Rabbit
Species Reactivity	Clostridium saccharobutylicum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium saccharobutylicum
Lead Time	Made-to-order (14-16 weeks)
Gene Names	fixB