



fliV Antibody

Product Code	CSB-PA334290XA01SBJ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P37588
Immunogen	Recombinant Salmonella muenchen fliV protein
Raised In	Rabbit
Species Reactivity	Salmonella muenchen
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Salmonella muenchen
Lead Time	Made-to-order (14-16 weeks)
Gene Names	fliV