



Tail spike Antibody

Product Code	CSB-PA312038XA01ECU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q04830
Immunogen	Recombinant Enterobacteria phage K1F (Bacteriophage K1F) Tail spike protein
Raised In	Rabbit
Species Reactivity	Enterobacteria phage K1F (Bacteriophage K1F)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Enterobacteria phage K1F (Bacteriophage K1F)
Lead Time	Made-to-order (14-16 weeks)
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.