



ORF125 Antibody

Product Code	CSB-PA137465XA01FQT
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7Y5D9
Immunogen	Recombinant Escherichia phage T5 (Enterobacteria phage T5) ORF125 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage T5 (Enterobacteria phage T5)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage T5 (Enterobacteria phage T5)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ORF125
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.