



ftt Antibody

Product Code	CSB-PA673208XA01MSI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q49610
Immunogen	Recombinant Methanopyrus kandleri (strain AV19 / DSM 6324 / JCM 9639 / NBRC 100938) ftt protein
Raised In	Rabbit
Species Reactivity	Methanopyrus kandleri (strain AV19 / DSM 6324 / JCM 9639 / NBRC 100938)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Methanopyrus kandleri (strain AV19 / DSM 6324 / JCM 9639 / NBRC 100938)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ftt
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.