



wbgU Antibody

Product Code	CSB-PA254821XA01FTW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7BJX9
Immunogen	Recombinant Plesiomonas shigelloides (Aeromonas shigelloides) wbgU protein
Raised In	Rabbit
Species Reactivity	Plesiomonas shigelloides (Aeromonas shigelloides)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Plesiomonas shigelloides (Aeromonas shigelloides)
Lead Time	Made-to-order (14-16 weeks)
Target Names	wbgU
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.