



gin Antibody

Product Code	CSB-PA654320XA01FSP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q38199
Immunogen	Recombinant Escherichia phage D108 (Bacteriophage D108) gin protein
Raised In	Rabbit
Species Reactivity	Escherichia phage D108 (Bacteriophage D108)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage D108 (Bacteriophage D108)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gin
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.