



gap Antibody

Product Code	CSB-PA009232XA01FUA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q59800
Immunogen	Recombinant Kitasatospora aureofaciens (Streptomyces aureofaciens) gap protein
Raised In	Rabbit
Species Reactivity	Kitasatospora aureofaciens (Streptomyces aureofaciens)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Kitasatospora aureofaciens (Streptomyces aureofaciens)
Lead Time	Made-to-order (14-16 weeks)
Target Names	gap
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.