



gyrB Antibody

Product Code	CSB-PA313952XA01MVW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P0C559
Immunogen	Recombinant Mycobacterium smegmatis gyrB protein
Raised In	Rabbit
Species Reactivity	Mycobacterium smegmatis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Mycobacterium smegmatis
Lead Time	Made-to-order (14-16 weeks)
Target Names	gyrB
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.