



gyrB Antibody

Product Code	CSB-PA340494XA01EYY
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P22840
Immunogen	Recombinant Proteus mirabilis gyrB protein
Raised In	Rabbit
Species Reactivity	Proteus mirabilis
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Proteus mirabilis
Lead Time	Made-to-order (14-16 weeks)
Gene Names	gyrB