



HAK1 Antibody

Product Code	CSB-PA343482XA01SDS
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P50505
Immunogen	Recombinant Schwanniomyces occidentalis (Yeast) (Debaryomyces occidentalis) HAK1 protein
Raised In	Rabbit
Species Reactivity	Schwanniomyces occidentalis (Yeast) (Debaryomyces occidentalis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Schwanniomyces occidentalis (Yeast) (Debaryomyces occidentalis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	HAK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.