



## grpE Antibody

|                     |                                                                                                 |
|---------------------|-------------------------------------------------------------------------------------------------|
| Product Code        | CSB-PA329828XA01BUD                                                                             |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                   |
| Uniprot No.         | P28609                                                                                          |
| Immunogen           | Recombinant Borrelia burgdorferi (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680) grpE protein |
| Raised In           | Rabbit                                                                                          |
| Species Reactivity  | Borrelia burgdorferi (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)                          |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                                    |
| Form                | Liquid                                                                                          |
| Conjugate           | Non-conjugated                                                                                  |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                |
| Purification Method | Antigen Affinity Purified                                                                       |
| Isotype             | IgG                                                                                             |
| Clonality           | Polyclonal                                                                                      |
| Product Type        | Polyclonal Antibody                                                                             |
| Immunogen Species   | Borrelia burgdorferi (strain ATCC 35210 / B31 / CIP 102532 / DSM 4680)                          |
| Lead Time           | Made-to-order (14-16 weeks)                                                                     |
| Target Names        | grpE                                                                                            |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                     |