



glmS Antibody

Product Code	CSB-PA902448XA01FTZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q05488
Immunogen	Recombinant Clostridium tetanomorphum glmS protein
Raised In	Rabbit
Species Reactivity	Clostridium tetanomorphum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium tetanomorphum
Lead Time	Made-to-order (14-16 weeks)
Gene Names	glmS