



licB Antibody

Product Code	CSB-PA767985XA01FUD
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q84C00
Immunogen	Recombinant Clostridium thermocellum (Ruminiclostridium thermocellum) licB protein
Raised In	Rabbit
Species Reactivity	Clostridium thermocellum (Ruminiclostridium thermocellum)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium thermocellum (Ruminiclostridium thermocellum)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	licB