



bglA Antibody

Product Code	CSB-PA338318XA01RMC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P45798
Immunogen	Recombinant Rhodothermus marinus (Rhodothermus obamensis) bglA protein
Raised In	Rabbit
Species Reactivity	Rhodothermus marinus (Rhodothermus obamensis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Rhodothermus marinus (Rhodothermus obamensis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	bglA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.