



celA Antibody

Product Code	CSB-PA345433XA01CLZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P54937
Immunogen	Recombinant Clostridium longisporum celA protein
Raised In	Rabbit
Species Reactivity	Clostridium longisporum
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Clostridium longisporum
Lead Time	Made-to-order (14-16 weeks)
Target Names	celA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.