



GSTY2 Antibody

Product Code	CSB-PA333503XA01FUF
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P30102
Immunogen	Recombinant Pichia kudriavzevii (Yeast) (Issatchenkia orientalis) GSTY2 protein
Raised In	Rabbit
Species Reactivity	Pichia kudriavzevii (Yeast) (Issatchenkia orientalis)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Pichia kudriavzevii (Yeast) (Issatchenkia orientalis)
Lead Time	Made-to-order (14-16 weeks)
Target Names	GSTY2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.