



hisJ Antibody

Product Code	CSB-PA152720XA01NFZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06758
Immunogen	Recombinant Neisseria gonorrhoeae hisJ protein
Raised In	Rabbit
Species Reactivity	Neisseria gonorrhoeae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Neisseria gonorrhoeae
Lead Time	Made-to-order (14-16 weeks)
Gene Names	hisJ