



hifA Antibody

Product Code	CSB-PA321214XA01HSZ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P14212
Immunogen	Recombinant Haemophilus influenzae hifA protein
Raised In	Rabbit
Species Reactivity	Haemophilus influenzae
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Haemophilus influenzae
Lead Time	Made-to-order (14-16 weeks)
Target Names	hifA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.