



13 Antibody

Product Code	CSB-PA293900XA01FTP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XQ99
Immunogen	Recombinant Escherichia phage T1 (Bacteriophage T1) 13 protein
Raised In	Rabbit
Species Reactivity	Escherichia phage T1 (Bacteriophage T1)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia phage T1 (Bacteriophage T1)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	13