



IplA Antibody

Product Code	CSB-PA487562XA01EOR
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	B7LNS6
Immunogen	Recombinant Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73) IplA protein
Raised In	Rabbit
Species Reactivity	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Escherichia fergusonii (strain ATCC 35469 / DSM 13698 / CDC 0568-73)
Lead Time	Made-to-order (14-16 weeks)
Target Names	IplA
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.