



LUXB2 Antibody

Product Code	CSB-PA341048XA01EUW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P29239
Immunogen	Recombinant Photobacterium leiognathi LUXB2 protein
Raised In	Rabbit
Species Reactivity	Photobacterium leiognathi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Species	Photobacterium leiognathi
Lead Time	Made-to-order
Gene Names	luxB