



lovD Antibody

Product Code	CSB-PA196527XA01DPC
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y7D1
Immunogen	Recombinant Aspergillus terreus lovD protein
Raised In	Rabbit
Species Reactivity	Aspergillus terreus
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aspergillus terreus
Lead Time	Made-to-order (14-16 weeks)
Gene Names	lovD