



luxD Antibody

Product Code	CSB-PA324166XA01EUW
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21309
Immunogen	Recombinant Photobacterium leiognathi luxD protein
Raised In	Rabbit
Species Reactivity	Photobacterium leiognathi
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Photobacterium leiognathi
Lead Time	Made-to-order (14-16 weeks)
Gene Names	luxD