



luxY Antibody

Product Code	CSB-PA324217XA01FTU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P21578
Immunogen	Recombinant Aliivibrio fischeri (Vibrio fischeri) luxY protein
Raised In	Rabbit
Species Reactivity	Aliivibrio fischeri (Vibrio fischeri)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Aliivibrio fischeri (Vibrio fischeri)
Lead Time	Made-to-order (14-16 weeks)
Target Names	luxY
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.