



ME1 Antibody

Product Code	CSB-PA013634XA01FUP
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P28227
Immunogen	Recombinant Anas platyrhynchos (Mallard) (Anas boschas) ME1 protein
Raised In	Rabbit
Species Reactivity	Anas platyrhynchos (Mallard) (Anas boschas)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Anas platyrhynchos (Mallard) (Anas boschas)
Lead Time	Made-to-order (14-16 weeks)
Target Names	ME1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.