




VACWR117 Antibody

Product Code	CSB-PA356237XA01VAI
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P04318
Immunogen	Recombinant Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR)) VACWR117 protein (49-128aa)
Species Reactivity	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Vaccinia virus (strain Western Reserve) (VACV) (Vaccinia virus (strain WR))
Target Names	VACWR117

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: VACWR117 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 36 kDa

Observed band size: 34 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.