



mdh Antibody

Product Code	CSB-PA013623XA01FUQ
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P80039
Immunogen	Recombinant Chlorobium tepidum (strain ATCC 49652 / DSM 12025 / NBRC 103806 / TLS) mdh protein
Raised In	Rabbit
Species Reactivity	Chlorobium tepidum (strain ATCC 49652 / DSM 12025 / NBRC 103806 / TLS)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Chlorobium tepidum (strain ATCC 49652 / DSM 12025 / NBRC 103806 / TLS)
Lead Time	Made-to-order (14-16 weeks)
Gene Names	mdh