



MD1 Antibody

Product Code	CSB-PA712456XA01NHE
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FSF0
Immunogen	Recombinant Nicotiana tabacum (Common tobacco) MD1 protein
Raised In	Rabbit
Species Reactivity	Nicotiana tabacum (Common tobacco)
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Nicotiana tabacum (Common tobacco)
Lead Time	Made-to-order (14-16 weeks)
Target Names	MD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.